



Fluo Sens Integrated – Fluorescence Detector



DIALUNOX GmbH, Bodenseeallee 20, D-78333 Stockach

ESM030-DF-1001

June 2021



1 Introduction

The Fluo Sens Integrated is highly suited for integration into compact and mobile fluorescence detection systems or for online measurements in automated processes.

Based on DIALUNOX' unique technology, the confocal optics allow contact free measurements on surfaces and in liquids, sometimes even in the presence of ambient light. The highly sensitive detector measures fluorescent samples that are present in low concentrations and is easy to integrate as a component into your instrument.

The compact and robust module contains a complete fluorescence measurement device comprising precise micro optics, powerful excitation light sources, highly sensitive sensors, no movable parts and microprocessor controlled electronics.

Calibration data and correction factors can be stored in the detector as software parameters and applied directly to the raw data. Measurement results can be accessed via a simple serial interface through external devices such as personal computers, embedded controllers (ECs) or programmable logic controllers (PLCs).

2 Technical specifications

Performance

Definition	Description	
Absolute lower detection limit LLOD	e.g. @ 470nm excitation and 520nm detection 15 pMol fluorescein-sodium in 0.1 M sodium hydroxide	
Dynamic range	3.6 orders of magnitude (0 2500 arb. u.* 1 linear range)	
	Adjustable range via Software by factor 10	
	Adjustable range via hardware by factor 1000	
Noise level	< 0.2% at max. range	
Excitation	High performance-LED with feedback loop for stabilization	
Detection	High amplification with high signal-to-noise ratio (between 98.5% and 99.95% depending on multiple factors), precision Si-Photodiode	
Measurement intervals	0.1 seconds to hours	
	Signal settling time ("On-Delay"): 300 ms	
	Signal fall time ("Off-Delay"): 300 ms	
	Measurement frequency in "scan mode": 200 Hz	
Detection area	1 mm ² to 25 mm ² (depends on working distance and used front lens)	
Distance (detector/object)	6 mm to 18 mm (depends on used front lens)	
Available excitation wavelengths**	365 nm 980 nm (two different spectral excitation ranges per detector)	
Available detection wavelengths**	460 nm 980 nm (two different spectral detection ranges per detector)	

* arb. u. = arbitrary unit

**Typical dyes: Alexa Fluor[®] 647, Pacific Blue EGFP, 5-FAM[™], ROX, HEX, Cy[®]3, Cy5, TAMRA, FITC



Environmental operating conditions

Definition	Description	
Temperature range	+10°C to +40°C	
Air humidity	20% - 70% rel. humidity, without condensation	
Air pressure	300 - 1060 hPa	

Mechanics

Definition	Description	
Housing	Aluminum	
Dimensions (without adapter)	64 x 47 x 17.8 mm	
Weight	< 90 g	

Electronics

Definition	Description	
Power supply	+5 V DC ±5%, ripple ≤ 20 mV	
Power input	Detector: 50 mA max	
	LED: ≤ 150 mA (depending on LED)	
Interface	Serial, 57600 baud, 1 start byte, 8 data bytes, no parity, 1 stop byte	
Interface models	TTL-level (3.3 V / 5 V tolerant)	
	RS232-level (±6 V)	
Connectors	10-way flex cable, MOLEX 98267-0257, adjustable length	

3 Physical dimensions





ESMO40-MB-xxxx detector dimensions (front side optical output)









 ${\tt ESMO30-MB-xxxx}\ detector\ dimensions\ (downside\ optical\ output).$

4 Connections





Overview

Pin number	Name	Description
1	Ground	Ground return for supply and communication signals
2	Ground	Ground return for supply and communication signals
3**	TxD	Output of the serial interface
4*	#Trigger	Low active input to trigger a measurement by hardware
5**	RxD	Input of the serial interface
6	Ground	Ground return for supply an communication signals
7*	#RESET	Low active input to reset the detector
		Detector supply voltage
8	+5 V	+5 V DC, ±5%, acceptable ripple <20 mV
		power consumption ≤40 mA
9	Ground	Ground return for supply an communication signals
10	VLED	LED supply voltage +5VDC

* These inputs do not provide internal pull up resistors. In normal operation they should be tied to +5 V/+3.3 V.

** The levels of the serial interface are with respect to the detector type either ±6V (RS232), 0/+3.3 V (TTL).

Visit www.fluorescence-measurement.com and discover more!

Not for use in diagnostic procedures.

For up-to-date licensing information and product-specific disclaimers, see the respective DIALUNOX kit handbook or user manual. DIALUNOX kit handbooks and user manuals are available at www.dialunox.com or can be requested from DIALUNOX Technical Services or your local distributor.

Trademarks: Alexa Fluor[®], FAM[™] (Life Technologies Corporation); Cy[®] (GE Healthcare). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

DIALUNOX GmbH / Bodenseeallee 20 / 78333 Stockach /Germany Ordering: sales@dialunox.com / Website: www.dialunox.com and www.fluorescence-measurement.com

© 2021 DIALUNOX, all rights reserved.